# **NAG Library Function Document**

# nag sparse sym basic diagnostic (f11gfc)

### 1 Purpose

nag\_sparse\_sym\_basic\_diagnostic (f11gfc) is the third in a suite of three functions for the iterative solution of a symmetric system of simultaneous linear equations (see Golub and Van Loan (1996)). nag\_sparse\_sym\_basic\_diagnostic (f11gfc) returns information about the computations during an iteration and/or after this has been completed. The first function of the suite, nag\_sparse\_sym\_basic\_set up (f11gdc), is a setup function, the second function, nag\_sparse\_sym\_basic\_solver (f11gec) is the proper iterative solver.

These three functions are suitable for the solution of large sparse symmetric systems of equations.

## 2 Specification

### 3 Description

nag\_sparse\_sym\_basic\_diagnostic (f11gfc) returns information about the solution process. It can be called both during a monitoring step of the solver nag\_sparse\_sym\_basic\_solver (f11gec), or after this solver has completed its tasks. Calling nag\_sparse\_sym\_basic\_diagnostic (f11gfc) at any other time will result in an error condition being raised.

For further information you should read the documentation for nag\_sparse\_sym\_basic\_setup (f11gdc) and nag sparse sym basic solver (f11gec).

#### 4 References

Golub G H and Van Loan C F (1996) *Matrix Computations* (3rd Edition) Johns Hopkins University Press, Baltimore

## 5 Arguments

1: **itn** – Integer \* Output

On exit: the number of iterations carried out by nag sparse sym basic solver (f11gec).

2: **stplhs** – double \* Output

On exit: the current value of the left-hand side of the termination criterion used by nag sparse sym basic solver (f11gec).

3: stprhs – double \* Output

On exit: the current value of the right-hand side of the termination criterion used by nag sparse sym\_basic\_solver (f11gec).

4: anorm – double \* Output

On exit: for CG and SYMMLQ methods, the norm  $||A||_1 = ||A||_{\infty}$  when either it has been supplied to nag\_sparse\_sym\_basic\_setup (f11gdc) or it has been estimated by nag\_sparse\_sym\_

Mark 26 f11gfc.1

f11gfc NAG Library Manual

basic\_solver (fl1gec) (see also Sections 3 and 5 in nag\_sparse\_sym\_basic\_setup (fl1gdc)). Otherwise, anorm = 0.0 is returned.

For MINRES method, an estimate of the infinity norm of the preconditioned matrix operator.

5: sigmax – double \*

Output

On exit: for CG and SYMMLQ methods, the current estimate of the largest singular value  $\sigma_1(\bar{A})$  of the preconditioned iteration matrix  $\bar{A} = E^{-1}AE^{-T}$ , when either it has been supplied to nag\_sparse\_sym\_basic\_setup (fl1gdc) or it has been estimated by nag\_sparse\_sym\_basic\_solver (fl1gec) (see also Sections 3 and 5 in nag\_sparse\_sym\_basic\_setup (fl1gdc)). Note that if its < itn then sigmax contains the final estimate. If, on final exit from nag\_sparse\_sym\_basic\_solver (fl1gec), its = itn, then the estimation of  $\sigma_1(\bar{A})$  may have not converged; in this case you should look at the value returned in sigerr. Otherwise, sigmax = 0.0 is returned.

For MINRES method, an estimate of the final transformed residual.

6: **its** – Integer \*

Output

On exit: for CG and SYMMLQ methods, the number of iterations employed so far in the computation of the estimate of  $\sigma_1(\bar{A})$ , the largest singular value of the preconditioned matrix  $\bar{A} = E^{-1}AE^{-T}$ , when  $\sigma_1(\bar{A})$  has been estimated by nag\_sparse\_sym\_basic\_solver (f11gec) using the bisection method (see also Sections 3, 5 and 9 in nag\_sparse\_sym\_basic\_setup (f11gdc)). Otherwise, its = 0 is returned.

7: **sigerr** – double \*

Output

On exit: for CG and SYMMLQ methods, if  $\sigma_1(\bar{A})$  has been estimated by nag\_sparse\_sym\_ba sic\_solver (f11gec) using bisection,

$$\mathbf{sigerr} = \max\left(\frac{\left|\sigma_1^{(k)} - \sigma_1^{(k-1)}\right|}{\sigma_1^{(k)}}, \frac{\left|\sigma_1^{(k)} - \sigma_1^{(k-2)}\right|}{\sigma_1^{(k)}}\right),$$

where k = its denotes the iteration number. The estimation has converged if  $sigerr \le sigtol$  where sigtol is an input argument to nag\_sparse\_sym\_basic\_setup (f11gdc). Otherwise, sigerr = 0.0 is returned.

For MINRES method, an estimate of the condition number of the preconditioned matrix.

8: work[lwork] - const double

Communication Array

On entry: the array work as returned by nag\_sparse\_sym\_basic\_solver (f11gec) (see also Section 3 in nag sparse sym basic solver (f11gec)).

9: **lwork** – Integer

Input

On entry: the dimension of the array work (see also Section 5 in nag\_sparse\_sym\_basic\_setup (f11gdc)).

Constraint: **lwork**  $\geq$  120.

**Note:** although the minimum value of **lwork** ensures the correct functioning of nag\_sparse\_sym\_basic\_diagnostic (f11gfc), a larger value is required by the iterative solver nag\_sparse\_sym\_basic\_solver (f11gec) (see also Section 5 in nag\_sparse\_sym\_basic\_setup (f11gdc)).

10: **fail** – NagError \*

Input/Output

The NAG error argument (see Section 2.7 in How to Use the NAG Library and its Documentation).

f11gfc.2 Mark 26

## 6 Error Indicators and Warnings

### NE\_ALLOC\_FAIL

Dynamic memory allocation failed.

See Section 2.3.1.2 in How to Use the NAG Library and its Documentation for further information.

#### NE\_BAD\_PARAM

On entry, argument  $\langle value \rangle$  had an illegal value.

### NE INT

On entry, **lwork** =  $\langle value \rangle$ . Constraint: **lwork** > 120.

### NE INTERNAL ERROR

An internal error has occurred in this function. Check the function call and any array sizes. If the call is correct then please contact NAG for assistance.

An unexpected error has been triggered by this function. Please contact NAG. See Section 2.7.6 in How to Use the NAG Library and its Documentation for further information.

### NE\_NO\_LICENCE

Your licence key may have expired or may not have been installed correctly. See Section 2.7.5 in How to Use the NAG Library and its Documentation for further information.

#### NE OUT OF SEQUENCE

nag sparse sym basic diagnostic (f11gfc) has been called out of sequence.

### 7 Accuracy

Not applicable.

### 8 Parallelism and Performance

nag\_sparse\_sym\_basic\_diagnostic (f11gfc) is not threaded in any implementation.

### 9 Further Comments

None.

## 10 Example

See Section 10 in nag\_sparse\_sym\_basic\_setup (fl1gdc).

Mark 26 fl1gfc.3 (last)